

**Avantor™ Performance Materials**

100 Matsonford Rd #200

Radnor, PA 19087 USA

Technical Service Department

Toll free 1 855 282 6867, select option 5

Outside U.S. 1 610 573 2600, select option 5

technical.service@avantorsciences.com

Date: September 15, 2023

Re: Freezing of Glacial Acetic Acid

**Physicochemical Properties:**

1. **Purity:** Each of Avantor Performance Materials, Inc. products listed as Glacial Acetic Acid has an assay, or purity, 99.5 % min. or better.

2. **Freezing / Melting:** The Freezing (Melting) Point (i.e., the temperature at which a pure substance has its solid and liquid phases in equilibrium at atmospheric pressure (1 atm)<sup>1</sup>) of Glacial Acetic Acid is 16.7 deg C (63 deg F)<sup>2,3</sup>. As such, exposure of any of the above Acetic Acid products to temperatures less than 16.7 deg C (63 deg F) results in *freezing*, or solidification, of the initially pure liquid product. Freezing of Glacial Acetic Acid is quite common during winter months, in climates having temperatures falling below these values.

3. **Performance:** Freezing of Glacial Acetic Acid is totally reversible within the reasonable working / storage temperature range on the condition that none of the initial liquid phase, the solidified product, or the 'reversed' liquid phase product was contaminated, especially by vapors of other organic chemicals. This is the case for an initially sealed bottle or container of Glacial Acetic Acid in its initially, pure liquid form.

This means that all product specifications, physicochemical properties, and performance characteristics of Glacial Acetic Acid products are totally maintained if it is initially liquid, then frozen (solidified), then melted (liquefied), within the temperature range of +40 to -20 deg C (+104 to -4 deg F).

4. **Recommendations for Melting** (Glacial Acetic Acid): Avantor has no single, formal 'set of instructions' or 'step-by-step procedure' for 'melting' any of the above Glacial Acetic Acid products which may have frozen. Nevertheless, the following procedures and recommendations should be followed:

**Important Note:** Be aware of, read about, and use, all appropriate safety and handling procedures, guidelines, and equipment indicated in the Avantor product MSDS, of Glacial Acetic Acid. Refer to additional standard chemical safety, handling, and applications references about Glacial Acetic Acid, if necessary.



**Avantor™ Performance Materials**  
Technical Service Department  
100 Matsonford Rd #200  
Radnor, PA 19087 USA

Method 1. Allow the frozen acetic acid, in its bottle or container, to thaw at room temperature (typically 20 - 25 deg C). The time for thawing may range from a few hours, to a few days, depending upon the extent to which the initial liquid phase has frozen.

Alternatively, to expedite the 'melting' or 'thawing' process of the Acetic Acid product,

Method 2. Subject the frozen acetic acid to 'controlled' heating (e.g., thermostat controlled heated water bath), to no higher than +35 °C (95° F). For safety, to prevent a hazardous build up of acetic acid vapor pressure, prior to heating, slightly loosen the cap of the Acetic Acid bottle or container. Perform the controlled heating procedure in a ventilated hood or ventilated laboratory area, away from chemicals incompatible with Acetic Acid, according to the product Material Safety Data Sheet (i.e., MSDS). Perform this procedure until a single, liquid phase of Acetic Acid, is observed in the bottle or container.

Following return of the frozen Acetic Acid product to a single, liquid phase, recommendation is to maintain the bottle or container at temperatures above 16.7 deg C (63 deg F), for short- or long- term storage.

If you have any questions or require additional information, please contact Technical Services.

---

**Prepared by the Technical Service Department**

While the above information is provided in good faith and believed to be accurate as of the date provided, Avantor makes no representations or warranties as to the accuracy or completeness of such information. All Avantor products are subject to Avantor's terms and conditions of sale including the limitations of liability contained therein and any contrary terms and conditions are expressly rejected. As Avantor has no control over purchasers' uses of its products, Avantor expressly disclaims all liability with respect to same.

***Revision History***

Rev. 0; Jan 25, 2023 – New Document

Rev. 1; September 8, 2023 – Updated to new template. Reviewed for correctness.