
Potassium Phosphate Dibasic

Product Regulatory Data Sheet

Section 1 – Product Information

Products Covered

<u>Brand</u>	<u>Product Code</u>	<u>Product Description</u>	<u>MOC* code</u>
J.T.Baker®	3222	Potassium Phosphate, Dibasic U.S.P., A.C.S.	R
J.T.Baker®	3250	Potassium Phosphate, Dibasic, U.S.P. Multi-Compendial	HR
J.T.Baker®	3251	Potassium Phosphate, Dibasic U.S.P.	R
J.T.Baker®	3255	Potassium Phosphate, Dibasic, Anhydrous Multi-Compendial	HR
J.T.Baker®	325A	Potassium Phosphate, Dibasic, U.S.P. Multi-Compendial	HR
J.T.Baker®	CH11	Potassium Phosphate, Dibasic, ChP Multi-Compendial	R
J.T.Baker®	3253	Potassium Phosphate, Dibasic, U.S.P. Multi-Compendial	R

*MOC = Management of Change

Section 2 – Manufacturing, Packaging and Release Site Information

The product codes 3250, 3253, 3255, and 325A, listed in Section 1, are manufactured under current Good Manufacturing Practices (cGMPs) as set forth by ICH Q7 and International Pharmaceutical Excipients Council (IPEC) guidelines.

The product codes 3222, 3251, and CH11, listed in Section 1, are manufactured according to current Good Manufacturing Practices (cGMPs) as set forth by International Pharmaceutical Excipients Council (IPEC) guidelines.

A number of the cGMP produced products that are sold by Avantor may not be originally manufactured at our sites. However, we perform the analytical and stability testing for these products and repackage the products where applicable. With ISO and cGMP procedures in place at our facilities, we can ensure, and take complete responsibility for, the traceability and quality of the finished, packaged product that we offer.

For J.T.Baker® and Macron Fine Chemicals™ brand products, the Original Manufacturer and address will be referenced on the Certificate of Analysis as an alpha or alpha-numeric **manufacturer code** rather than listing the full name and address. This practice is compliant with both ICH Q7 Good Manufacturing Guidance for Active Pharmaceutical Ingredients (APIs) and IPEC guidelines and it meets

cGMP requirements. For instructions to decipher the manufacturer reference code please consult the Avantor website. Instructions can be found by visiting the Ask Avantor link under the Resources tab or by directly linking to www.askavantor.com Keyword: Manufacturer Code. Additional information on Avantor suppliers may be available under NDA. Please reach out to the support contact in Section 7 for additional supplier information inquiries.

Section 3 – Physical/Chemical Information

CAS #: 7758-11-4

Manufacturing Process: Synthesis- Chemical Reaction. Additional manufacturing process information may be disclosed under NDA upon request from the support contact in Section 7.

Raw Material Origin: Chemical

Section 4 – Regulatory Information

DMF: Avantor may hold Master File(s) for specified product codes, dependent on the country of interest. Inquire with the support contact in Section 7 for additional details.

BSE/TSE Status: The subject materials are manufactured from raw materials that contain NO animal parts, products, and/or by-products nor do they come in contact with animal parts, products, and/or by-products.

Allergen/Hypersensitivities Information: To the best of our knowledge, the allergens listed in the [US FDA](#), [EU Directive 2003/89/EC](#), and [TGO-91/92](#) are not known additives, by products, intermediate parts, or otherwise intentionally added during the manufacturing processes of the product.

Avantor does not produce any of the following types of products: antibiotics, penicillin, semi-synthetic penicillins, cephalosporins, other beta-lactams, cytotoxics, steroids, medicated feeds, or pesticides.

This product is manufactured using cGMP guidelines which provide controls that allow no potential for cross contamination of any allergens or other contaminants including aflatoxins. However, this product is not tested for the presence of these or any other allergens by Avantor, therefore, we do not have confirmation for the absence of any allergens in the product.

GMO Information: The subject materials, including any raw materials and processing aids, are NOT subject to genetic modification.

Residual Solvents/Organic Volatile Impurities (OVI) Information: The subject materials (all lots) comply with the requirements of the ICH Q3C Residual Solvents Guideline and USP <467> Residual

Solvents. No Class 1, 2, 3 or other solvents are used or produced in the manufacturing or purification of the product.

Elemental Impurities: Please see attached summary for Elemental Impurity information for listed products.

Kosher Status: For J.T.Baker® and Macron Fine Chemicals™ brand products, kosher certification is aligned to the Avantor packaging site as indicated on the product Certificate of Analysis. Please refer to the site-specific kosher certificate available on AskAvantor for our most up to date listing of kosher products at (www.askavantor.com Keyword: kosher).

Halal Status: For J.T.Baker® and Macron Fine Chemicals™ brand products, halal certification is aligned to the Avantor packaging site as indicated on the product Certificate of Analysis. Please refer to the site-specific halal certificate available on AskAvantor for our most up to date listing of halal products at (www.askavantor.com Keyword: halal).

GRAS Status: The United States Food and Drug Administration (FDA) have acknowledged that some chemicals may be considered Substances Generally Recognized as Safe (GRAS) in foods when used in accordance with the requirements and limitations per specific 21 CFR regnums. For the latest information on whether or not an Avantor product is considered GRAS, please visit the [Electronic Code of Federal Regulations](#).

Section 5 – Miscellaneous Product Information

Certificate of Analysis Date Format: The Manufactured Date and Expiration/Retest Date on the Certificate of Analysis are reported as YYYY-MM-DD. For example, the Manufactured Date for October 1, 2021 would be reported as 2021-10-01.

Lot Numbering System and Batch Description: For J.T.Baker® and Macron Fine Chemicals™ brand products, please refer to Ask Avantor for information concerning our lot/batch numbering system. (www.askavantor.com Keyword: Lot Number).

Batch Definition: A "batch" is a homogeneous unit of production; each batch of is from one single batch of the source supplier.

Shelf-Life Information: If a product has an assigned expiration or retest period, the date will appear on the Certificate of Analysis. For products that do not have assigned dates, please reach out to the support contact in Section 7 for additional stability inquiries.

Management of Change: For J.T.Baker® and Macron Fine Chemicals™ brand products, please refer to Management of Change link under the Working with Avantor tab on the Avantor website.

Country of Origin Statement: Country of Origin is indicated on the product Certificate of Analysis. If you require further documentation, please reach out to the Trade Compliance support contact in Section 7.

Storage Requirements: Please refer to the product's Certificate of Analysis or Product Specifications. In the absence of specific storage conditions listed on its specification sheet or Certificate of Analysis, products are to be stored in ambient conditions of temperature and humidity. We do not formally tie any specific temperature or humidity range with the "ambient" storage designation, but an example of a common temperature interpretation is 15-30°C. Our products are also packaged to protect from the normal variation in humidity during storage and shipment. Further handling and storage information may be found in Section 7 of the product's SDS sheet.

Certificates of Analysis: For J.T.Baker® and Macron Fine Chemicals™ brand products, please see the current list of product specifications using the Certificate/SDS Search tool on our website [here](#).

Safety Data Sheet: For J.T.Baker® and Macron Fine Chemicals™ brand products, please see the current product safety information using the Certificate/SDS Search tool on our website [here](#).

Avantor Site Certifications: Please see the current Avantor site certifications on our website [here](#).

Site Quality Overview: Avantor maintains a self-assessment modeled after IPEC guidelines which describes site and quality system information to support the manufacturing activities of this product. Please reach out to the support contact in Section 7 for a current copy of the Site Quality Overview.

Packaging Information: Please reach out to the support contact in Section 7 for current packaging specifications.

Section 6 – Revision History

Rev. 0; Oct. 1, 2007 – IPEC EIP format

Rev. 1; Oct. 22, 2008 – Section 4: updated residual solvents information

Rev. 2; Jul. 22, 2009 – Section 1: added product number 7787; Section 7: reviewed and updated contact information; Entire Document: new letterhead, added 'keywords' to Q&A Center of the customer support section of our web site references. (KES)

Rev. 3; Jul. 22, 2011 – Entire document: new letterhead and changed all references of "AskMBI" to "AskAvantor." Updated website links for new website; minor formatting; Section 1: changed Mallinckrodt to Macron brand; added MOC codes; Section 2: added GMP statement; Section 4: expanded Allergens list; added GRAS statement; Section 7: updated contact information. (PH, JLW)

Rev. 4; June 24, 2013 – Entire document: updated headquarters address; Section 4: added add'l allergens as listed in EU Directive 2003/89/EC; updated Residual Metallic Catalysts statement; separated Kosher/Halal status and added certification statement Section 5: added Management of Change information; Added COA Date Format statement; Section 7: removed contact list table and added CS/TS contact information(MCH)

Rev. 5- March 13, 2014- Updated EMEA statement to reflect current guideline revision. (MCH)

Rev. 6- Dec. 4, 2014 – Section 1: Added code 3222. (MCH)

Rev. 7; November 14, 2018 – Entire document: new letterhead (logo & headquarters address, updated email from @avantormaterials.com to @avantorinc.com. Section 1: Removed product code 7787. updated MOC code from R to HR for 3250; Section 4: Updated DMF statement Replaced Residual Metallic Catalyst statement with Elemental Impurities statement. Section 5: added Storage Requirement, Batch Definition, and Country of Origin Statements (MCH)

Rev. 8; April 5, 2021 – Entire Document: Minor formatting. Changed website and email address from avantorinc.com to avantorsciences.com; Section 1: Added product codes 3255 per NSPU-2644 and 325A per NPSU-2797; Section 2: Updated language to specify guideline used for the manufacture of each product code; Section 4: Updated DMF and Allergen/Hypersensitivities Information statements Updated Elemental Impurity report to include new product codes; Section 5: Removed Nutritional/Supplement Facts Labeling and Organic Status statements. Products are not F.C.C. (KH)

Rev. 9; June 28, 2023- Entire Document: Formatting. Section 1: Added product codes CH11 in accordance with BOM-008233 and 3253; Section 2: Minor formatting to verbiage; Section 3: Minor formatting to verbiage; Section 4: Remove Compendial Compliance statement. Minor formatting to DMF, Allergen/Hypersensitivities Information, Kosher Status, Halal Status, and GRAS Status statement verbiage; Section 5: Minor formatting to Certificate of Analysis Date Format, Shelf-Life Information, Management of Change, Country of Origin, Storage Requirements statements verbiage. Added statements for Certificates of Analysis, Safety Data Sheet, Avantor Site Certifications, Site Quality Overview, Packaging Information; Section 7: Updated Contact Information; Disclaimer: Minor formatting of verbiage. Updated Elemental Impurity report to include new product codes. (KH)

Rev. 10; July 11, 2023- Corrected typographical error for MOC code for 3253 from HR to R. (KH)

This electronic document is valid without a signature.

Section 7 – Contact Information

Technical Service

Phone: 1-855-282-6867 and 1-610-573-2600 (outside U.S.), select option 5

Email: Technical.Service@avantorsciences.com

Regulatory Support

Email: regulatory.support@avantorsciences.com

Trade Compliance

Email: Trade.Compliance@avantorsciences.com

While the above information is provided in good faith and believed to be accurate as of the date provided, Avantor makes no representations or warranties as to the accuracy or completeness of such information. All Avantor products are subject to Avantor's terms and conditions of sale including the limitations of liability contained therein and any contrary terms and conditions are expressly rejected. As Avantor has no control over purchasers' uses of its products, Avantor expressly disclaims all liability with respect to same.

The most current revision of this document is maintained on our website. Reviews and revisions are performed as warranted due to product changes or as part of the supplier audit cycle and managed under a validated document control system.

Avantor, Inc.
100 Matsonford Rd., Suite 200
Radnor, PA 19087 USA
www.avantorsciences.com

Material Name: Potassium Phosphate Dibasic **Product codes:** 3250, 3222, 3251, 3255, 325A 3253 CH11 **Date:** September 28, 2016 Rev. 3

Source/Type of Excipient: ☐ Mineral; ☐ Mineral derived; ☐ Plant; ☐ Plant derived; ☒ Synthetic; ☐ Fermentation derived

Other (explain):

Elemental Impurity		Class	Likely to be Present			If Known, Please Identify the Expected Concentration /Units (or Range)	Analytical Method Used (and Limit of Detection if Available)	Comments regarding source of information (i.e.; number of lots tested, frequency of testing, process understanding, etc.)
Arsenic (inorganic)	As	1	Yes <input type="checkbox"/>	No <input checked="" type="checkbox"/>	Unknown <input type="checkbox"/>	<0.05 ppm	ICP-MS (MRL=0.05 ppm)	Avg. of 3 batches
Cadmium	Cd	1	Yes <input type="checkbox"/>	No <input checked="" type="checkbox"/>	Unknown <input type="checkbox"/>	<0.05 ppm	ICP-MS (MRL=0.05 ppm)	Avg. of 3 batches
Mercury (inorganic)	Hg	1	Yes <input type="checkbox"/>	No <input checked="" type="checkbox"/>	Unknown <input type="checkbox"/>	<0.05 ppm	ICP-MS (MRL=0.05 ppm)	Avg. of 3 batches
Lead	Pb	1	Yes <input type="checkbox"/>	No <input checked="" type="checkbox"/>	Unknown <input type="checkbox"/>	<0.05 ppm	ICP-MS (MRL=0.05 ppm)	Avg. of 3 batches
Cobalt	Co	2A	Yes <input type="checkbox"/>	No <input checked="" type="checkbox"/>	Unknown <input type="checkbox"/>	<0.05 ppm	ICP-MS (MRL=0.05 ppm)	Avg. of 3 batches
Nickel	Ni	2A	Yes <input checked="" type="checkbox"/>	No <input checked="" type="checkbox"/>	Unknown <input type="checkbox"/>	0.5 ppm	ICP-MS (MRL=0.05 ppm)	Avg. of 3 batches
Vanadium	V	2A	Yes <input type="checkbox"/>	No <input checked="" type="checkbox"/>	Unknown <input type="checkbox"/>	<0.05 ppm	ICP-MS (MRL=0.05 ppm)	Avg. of 3 batches
Silver	Ag	2B	Yes <input type="checkbox"/>	No <input checked="" type="checkbox"/>	Unknown <input type="checkbox"/>	<0.05 ppm	ICP-MS (MRL=0.05 ppm)	Avg. of 3 batches
Gold	Au	2B	Yes <input type="checkbox"/>	No <input checked="" type="checkbox"/>	Unknown <input type="checkbox"/>	<0.05 ppm	ICP-MS (MRL=0.05 ppm)	Avg. of 3 batches

Elemental Impurity		Class	Likely to be Present			If Known, Please Identify the Expected Concentration /Units (or Range)	Analytical Method Used (and Limit of Detection if Available)	Comments regarding source of information (i.e.; number of lots tested, frequency of testing, process understanding, etc.)
			Yes <input type="checkbox"/>	No <input checked="" type="checkbox"/>	Unknown <input type="checkbox"/>			
Iridium	Ir	2B	Yes <input type="checkbox"/>	No <input checked="" type="checkbox"/>	Unknown <input type="checkbox"/>	<0.05 ppm	ICP-MS (MRL=0.05 ppm)	Avg. of 3 batches
Osmium	Os	2B	Yes <input type="checkbox"/>	No <input checked="" type="checkbox"/>	Unknown <input type="checkbox"/>	<0.05 ppm	ICP-MS (MRL=0.05 ppm)	Avg. of 3 batches
Palladium	Pd	2B	Yes <input type="checkbox"/>	No <input checked="" type="checkbox"/>	Unknown <input type="checkbox"/>	<0.05 ppm	ICP-MS (MRL=0.05 ppm)	Avg. of 3 batches
Platinum	Pt	2B	Yes <input type="checkbox"/>	No <input checked="" type="checkbox"/>	Unknown <input type="checkbox"/>	<0.05 ppm	ICP-MS (MRL=0.05 ppm)	Avg. of 3 batches
Rhodium	Rh	2B	Yes <input type="checkbox"/>	No <input checked="" type="checkbox"/>	Unknown <input type="checkbox"/>	<0.05 ppm	ICP-MS (MRL=0.05 ppm)	Avg. of 3 batches
Ruthenium	Ru	2B	Yes <input type="checkbox"/>	No <input checked="" type="checkbox"/>	Unknown <input type="checkbox"/>	<0.05 ppm	ICP-MS (MRL=0.05 ppm)	Avg. of 3 batches
Selenium	Se	2B	Yes <input type="checkbox"/>	No <input checked="" type="checkbox"/>	Unknown <input type="checkbox"/>	<0.05 ppm	ICP-MS (MRL=0.05 ppm)	Avg. of 3 batches
Thallium	Tl	2B	Yes <input type="checkbox"/>	No <input checked="" type="checkbox"/>	Unknown <input type="checkbox"/>	<0.05 ppm	ICP-MS (MRL=0.05 ppm)	Avg. of 3 batches
Barium	Ba	3	Yes <input type="checkbox"/>	No <input checked="" type="checkbox"/>	Unknown <input type="checkbox"/>	<0.05 ppm	ICP-MS (MRL=0.05 ppm)	Avg. of 3 batches
Chromium	Cr	3	Yes <input checked="" type="checkbox"/>	No <input type="checkbox"/>	Unknown <input type="checkbox"/>	0.6 ppm	ICP-MS (MRL=0.05 ppm)	Avg. of 3 batches
Copper	Cu	3	Yes <input type="checkbox"/>	No <input checked="" type="checkbox"/>	Unknown <input type="checkbox"/>	<0.05 ppm	ICP-MS (MRL=0.05 ppm)	Avg. of 3 batches
Lithium	Li	3	Yes <input type="checkbox"/>	No <input checked="" type="checkbox"/>	Unknown <input type="checkbox"/>	<0.05 ppm	ICP-MS (MRL=0.05 ppm)	Avg. of 3 batches
Molybdenum	Mo	3	Yes <input type="checkbox"/>	No <input checked="" type="checkbox"/>	Unknown <input type="checkbox"/>	<0.05 ppm	ICP-MS (MRL=0.05 ppm)	Avg. of 3 batches

Elemental Impurity		Class	Likely to be Present			If Known, Please Identify the Expected Concentration /Units (or Range)	Analytical Method Used (and Limit of Detection if Available)	Comments regarding source of information (i.e.; number of lots tested, frequency of testing, process understanding, etc.)
			Yes <input checked="" type="checkbox"/>	No <input type="checkbox"/>	Unknown <input type="checkbox"/>			
Antimony	Sb	3	Yes <input checked="" type="checkbox"/>	No <input type="checkbox"/>	Unknown <input type="checkbox"/>	0.14 ppm	ICP-MS (MRL=0.05 ppm)	Avg. of 3 batches
Tin	Sn	3	Yes <input type="checkbox"/>	No <input checked="" type="checkbox"/>	Unknown <input type="checkbox"/>	<0.05 ppm	ICP-MS (MRL=0.05 ppm)	Avg. of 3 batches

Reference: ICH Q3D Guideline for Elemental Impurities, Step 4 version, September 2014



David L. Cugini, Sr. QA Analyst

Prepared by the Technical Service Department

Avantor™ Performance Materials, LLC

Avantor™ Performance Materials, LLC provides the information contained herein in good faith but makes no representations or warranties, either expressed or implied, including without limitation any warranties of merchantability or fitness for a particular purpose with respect to the information set forth herein or the product(s) to which the information refers. Accordingly, Avantor™ Performance Materials, LLC will not be responsible for damages resulting from the use of or reliance upon this information.

Trademarks are owned by Avantor Performance Materials, LLC or its affiliates unless otherwise noted.